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[EN] Consumer Information

NASCAR is a hypoallergenic, transparent, rapidly drying silicone gel for the treatment of scars of different origin.

Gel NASCAR is easily applied to all areas of the skin, including exposed areas like the face and neck; as well as joints and hairy areas without the need for shaving.

When used as directed gel rapidly dries to form a silicone gel sheet. A protective layer, that is gas permeable and waterproof, hydrates and protects a scar from chemical and microbial invasion, allowing scar to normalize the collagen synthesis cycle.

Gel NASCAR is used to soften and flatten scars, relieve the itching and discomfort of scars, redness and discoloration associated with scars. Gel NASCAR is used for treatment of all types of scars including those resulting from: general surgery, trauma, chronic wounds, cesarean operation, burns, bites, acne and other diseases of the skin.

In the International Clinical Recommendations on scar management stated that «silicone gel should be used as first line therapy in the initial management of scars and in the prevention of hypertrophic and keloid scars».*

Indications for use

Gel NASCAR is used for the prevention and treatment of all types of scars, including keloid and hypertrophic scars. NASCAR could be used in combination with other invasive methods for improving of overall result, including surgical aggression and sealing, injections and pressure garments. Gel NASCAR is suitable for skin of children of all ages, for women during pregnancy and lactation; for sensitive and allergyprone skin.

In order to ensure best results it is recommended to start the use of NASCAR as soon as possible after the period of wound epithelialization. After surgical measure with stitching NASCAR is applied from 10 day after surgery.

Directions for use

Ensure that the scar is clean and dry. Apply a thin layer of gel NASCAR to the scar and allow to dry. Gel should be applied once daily, or twice daily to exposed areas or as advised by a physician. For best results gel should be maintained with in continuous contact with the skin. When applied correctly, gel should be dry in 5–7 minutes. If it takes longer to dry you have probably applied too much. Gently remove the excess with a tissue and allow the drying process of gel to continue. Once dry, gel can be covered by sunscreen, cosmetics, pressure garments or casts.

Recommended duration of treatment

NASCAR is recommended for use during at least 60–90 days (24/7). Larger and older scars may take longer. Continued use is recommended until no further improvement is seen.

How much gel NASCAR do I need?

- 5 g of gel NASCAR is enough to treat a 2–3 cm scar for at least 60–90 days.
- 15 g of gel NASCAR is enough to treat a 7–10 cm scar for at least 60–90 days.
- 50 g of gel NASCAR is enough to treat an area of 6×12 cm² for at least 60–90 days.

Warnings

- Gel should not be applied on open wounds.
- Gel should not be placed in contact with eyes or mucous membranes.
- Gel should not be applied over other skin treatments without the advice of a physician.
- Gel may stain clothing if not completely dry.

For external use only. If irritation occurs, discontinue

use and consult physician. Keep away from sources of heat in the original package, out of the reach of children. Store at temperature below 25 °C.

Do not use after the expiration date, which printed on the package.

Composition

1 g of gel contains: 400 mg of dimethicone, 400 mg of dimethicone and dimethiconol, 200 mg of dimethicone and dimethicone crosspolymer.

Package

15 g or 50 g in tube. 1 tube in carton box.

15 g in roll-on bottle. 1 bottle in carton box.

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Amaha Pharma LTD, 72 Hammersmith Road, London, W14 8TH, United Kingdom